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LIGHTED SUCTION FOR CERUMEN REMOVAL - USER INSTRUCTIONS

Assembly/Disassembly:

1. To assemble the Lighted Suction Device, insert the blunt end of the Lighted Suction Handle into the light source by aligning the pegs with the grooves inside the light source. Push down and twist **FIRMLY** clockwise until rotation stops. This will lock the handle into place and activate the light.
2. Attach the suction hose (extending from pump) to the barbed hose adapter of the handle by pushing firmly.
3. To attach the magnification lens, snap the lens onto the neck of the blue light source so that the top of the lens angles toward the Lighted Suction Handle.
4. To disassemble, push the Lighted Suction Handle into the light source and rotate counter clockwise to remove. The light will go off when the tip is removed.

The Lighted Suction Handle is intended for single-patient use. The light source and magnification lens are designed for multi-procedure use. A light source and lens are provided in each box of product.

Instructions for Use:

1. Examine the ear using an otoscope noting the location and depth of cerumen or foreign body.
2. Under the direct visualization of the ear canal (using the Lighted Suction Handle with or without magnification) carefully approach the cerumen with your thumb removed from the fine tuning suction port on the handle.
3. Once in contact with the cerumen, adjust suction level by gradually covering the opening of the suction port until cerumen is removed or grasped against the suction tip. **NOTE:** Keep a cup of hot water nearby in case the cerumen suction handle becomes clogged. Drawing hot water into the Lighted Suction Handle may help clear the clog.
4. Once the procedure is complete, remove the Lighted Suction Handle from the light source and tubing then dispose of handle. Save the magnification lens and light source for your next procedure.

 **CAUTION:** Discontinue cerumen/foreign body removal immediately if bleeding, irritation, or other trauma to the ear canal or tympanic membrane occurs.

 **CAUTION:** This product is a spring-loaded device; improper assembly may cause ejection of the Lighted Suction Handle.

 **CAUTION:** DO NOT hold device by light source as accidental ejection of the Lighted Suction Handle may occur if not assembled properly.

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©2024 Bionix LLC | 1670 Indian Wood Circle