Disposable Nasal Speculum

INSTRUCTIONS FOR USE

R_x ONLY

MADE IN THE USA

**Importer, Importeur, Importör,
Importador, Importateur,
Importeur, Importør, Dovoze,
Maahantuoja, Importatore,
Εισαγωγέας, Importuotojas:**

Obelis s.a
Boulevard Général Wahis 53
1030 Brussels, BELGIUM
Tel: +(32) 2.732.59.54
Fax: +(32) 2.732.60.03
E-Mail: mail@obelis.net

Bionix LLC
1670 Indian Wood Circle
Maumee, Ohio 43537, United States
Phone: +1 419.727.8421
Fax: +1 419.727.8426
www.Bionix.com

MedEnvoy

Prinses Margrietplantsoen 33 - Suite 123
2595 AM The Hague
The Netherlands

©2026 Bionix LLC | 1670 Indian Wood Circle | Maumee, Ohio 43537 | United States
Phone: 800.551.7096 (+1 419.727.8421) | Fax: 800.455.5678 (+1 419.727.8426)
Email: bionix@bionix.com (international@bionix.com) | www.bionix.com
RM95-0028 | 2026/07/10 | Rev. 02

Disposable Nasal Speculum Instructions for Use

English	3
Deutsch (German)	4
Español (Spanish)	5
Français (French)	6
Svenska (Swedish)	7
Polski (Polish)	8
Česky (Czech)	9
Português (Portuguese)	10
ελληνικά (Greek)	11
Nederlands (Dutch)	12
Suomi (Finnish)	13
Italiano (Italian)	14
Lietuvių k. (Lithuanian)	15
Dansk (Danish)	16
Bokmål (Norwegian)	17
Română (Romanian)	18

ENGLISH

DISPOSABLE NASAL SPECULUM - USER INSTRUCTIONS

Intended Use: The intended use of the Disposable Nasal Speculum is to allow direct examination of the anterior part of the nose.

Indications for Use: The Disposable Nasal Speculum is indicated for use in diagnosis, extraction of foreign bodies, and nasal packing by healthcare professionals on all patient populations.

Instructions:

To use, hold the Disposable Nasal Speculum comfortably in one hand. Position the patient in a manner that ensures proper access and adequate lighting. Gently insert the Disposable Nasal Speculum tip into the patient's nares, with the blades in the "closed" position. Slowly squeeze the handle to spread the Disposable Nasal Speculum blades and open the nares. Stop when adequate exposure has been obtained.

⚠ CAUTION: Discontinue use of the Disposable Nasal Speculum if the patient experiences pain or trauma to the nose or the nasal mucosa.

⚠ WARNING: Cross-contamination risk. Do not reuse the Disposable Nasal Speculum as this may spread contamination from one patient to another patient.

Medical Device Reporting: Notice to Users and/or Patients in EU: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Notice for Users and/or Patients in Switzerland: Any serious incident that has occurred in relation to the device should be reported without delay to Bionix at bionix@bionix.com and to Swissmedic (www.swissmedic.ch).

Contraindications: There are no absolute contraindications related to the use of the Disposable Nasal Speculum.

Residual Risk: Risk associated with the use of this product has been reduced as far as possible, but the product cannot completely eliminate potential patient or user harm arising from the following:

- Harm from mechanical hazards
- Harm from misuse, or use error
- Harm from unanticipated origins

	Medical Device		Do Not Use if Package is Damaged
	Manufacturer		Do Not Reuse Single-Use Device
	Date of manufacture		Prescription Only or "For Use by or on the order of a licensed medical professional"
	Lot Code		Consult Instructions For Use (IFU)
	Reorder Number		Caution
	European Conformity		Warning
	Authorized representative in the European Community		

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