ShotBlocker®

INSTRUCTIONS FOR USE

Obelis s.a
Boulevard Général Wahis 53
1030 Brussels, BELGIUM
Tel: +(32) 2.732.59.54
Fax: +(32) 2.732.60.03
E-Mail: mail@obelis.net

R_x ONLY

Bionix LLC
1670 Indian Wood Circle
Maumee, Ohio 43537, United States
Phone: +1 419.727.8421
Fax: +1 419.727.8426
www.Bionix.com

MADE IN THE USA

**Importer, Importeur, Importör,
Importador, Importateur,
Importeur, Importør, Dovožce,
Maahantuoja, Importatore,
Εισαγωγέας, Importuotojas:**

MedEnvoy

Prinses Margrietplantsoen 33 - Suite 123
2595 AM The Hague
The Netherlands

©2026 Bionix LLC | 1670 Indian Wood Circle | Maumee, Ohio 43537 | United States
Phone: 800.551.7096 (+1 419.727.8421) | Fax: 800.455.5678 (+1 419.727.8426)
Email: bionix@bionix.com (international@bionix.com) | www.bionix.com
RM95-0031 | 2026/07/10 | Rev. 03

ShotBlocker® Instructions for Use

English	3
Deutsch (German)	4
Español (Spanish)	5
Français (French)	6
Svenska (Swedish)	7
Polski (Polish)	8
Česky (Czech)	9
Português (Portuguese)	10
ελληνικά (Greek)	11
Nederlands (Dutch)	12
Suomi (Finnish)	13
Italiano (Italian)	14
Lietuvių k. (Lithuanian)	15
Dansk (Danish)	16
Bokmål (Norwegian)	17
Română (Romanian)	18
Magyar (Hungarian)	19

ENGLISH

SHOTBLOCKER® - USER INSTRUCTIONS

Intended use: ShotBlocker is intended to distract, comfort, and reduce the anxiety of a patient receiving minor injections and/or needle pokes.

Indications for use: ShotBlocker is designed to distract and comfort patients receiving intra-muscular or subcutaneous injections and needle pokes.

Instructions:

1. Select the injection site and prep the skin as usual.
2. Hold ShotBlocker so that the blunt contact points touch the patient's skin at the injection site.
3. Press ShotBlocker **FIRMLY** against the skin.
Note: DO NOT MOVE OR REMOVE SHOTBLOCKER UNTIL THE INJECTION HAS BEEN COMPLETED.
4. Immediately administer the injection in the usual manner through or near the central opening of ShotBlocker. For subcutaneous injections, angle the needle as needed to give the injection. Note: IF MORE THAN 20 SECONDS ELAPSE BETWEEN THE PLACEMENT OF SHOTBLOCKER AND THE INJECTION, COMPLETELY REMOVE SHOTBLOCKER FROM THE SKIN. REPEAT THE PROCESS BEGINNING WITH STEP 2.

Helpful Hints When Using ShotBlocker®

- Explain how the ShotBlocker will use the distraction technique to hide the pain.
- Calm patients by letting them feel the device and demonstrate how it works.
- Have the patient look away while you perform the injection.
- Press the ShotBlocker firmly onto the patient making sure the contact points touch the skin. (This creates the necessary sensation for ShotBlocker to be effective)
- Best results are achieved when performed on children ages 4 and up.
- ShotBlocker significantly reduces the pain and discomfort associated with routine injections 80% of the time. (Maumee Pediatric Associates, 2002)
- The ShotBlocker does not eliminate the sting from medications but distracts from the pain of the needle poke.

Disposal

ShotBlocker should be recycled or disposed according to current local, regional, state, and/or federal laws and regulations pertaining to the safe disposal of medical devices and accessories. If the user is uncertain how to proceed, they should contact Bionix, LLC for guidance on safe disposal protocols.

WARNING: Cross-contamination risk.

The ShotBlocker is a single-use device – do not reuse as this may spread contamination from one patient to another patient.

Medical Device Reporting: Notice to Users and/or Patients in EU: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Notice for Users and/or Patients in Switzerland:

Any serious incident that has occurred in relation to the device should be reported without delay to Bionix at bionix@bionix.com and to Swissmedic (www.swissmedic.ch).

Contraindications: There are no absolute contraindications for use of the ShotBlocker.

Residual Risk: Risk associated with the use of this product has been reduced as far as possible, but the product cannot completely eliminate potential patient or user harm arising from the following:

- Harm from mechanical hazards
- Harm from misuse, or use error
- Harm from unanticipated origins

	Medical Device		Do Not Use if Package is Damaged
	Manufacturer		Do Not Reuse Single-Use Device
	Date of manufacture		Prescription Only or "For Use by or on the order of a licensed medical professional"
	Lot Code		Consult Instructions For Use (IFU)
	Reorder Number		Caution
	European Conformity		Warning
	Authorized representative in the European Community		

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