Igloo® Wound Irrigation Shield

INSTRUCTIONS FOR USE

STERILE	R
---------	---

R_x ONLY

MADE IN THE USA

©2026 Bionix LLC | 1670 Indian Wood Circle | Maumee, Ohio 43537 | United States
Phone: 800.551.7096 (+1 419.727.8421) | Fax: 800.455.5678 (+1 419.727.8426)
Email: bionix@bionix.com (international@bionix.com) | www.bionix.com
RM95-0505 | 2026/07/10 | Rev. 01

Igloo® Wound Irrigation Shield IFU

English	3
Français (French)	4

ENGLISH

Intended use: The Igloo Wound Irrigation Shield is intended to be used during wound irrigation to protect the user from backslash as well as controlling the effluent runoff.

Indications for use: The Igloo Wound Irrigation Shield is indicated for use during the general cleansing of wound sites. The Igloo Wound Irrigation Shield may be used on **all ages of patients**.

USER INSTRUCTIONS:

1. Prep and drape the wound in a sterile fashion. Prepare a basin with sterile saline or other desired irrigant solution.
2. Fill a standard luer lock syringe or other delivery device with the irrigation solution.
3. Attach the Igloo Wound Irrigator to the tip of the syringe or other delivery device.
4. Place the Igloo Wound Irrigator over the wound, lightly touching the skin of the patient. For large, spacious (or chronic) wounds, place Igloo just above the wound surface.
5. Irrigate the wound in the usual fashion; syringe may need to be refilled one or more times as needed.
6. After use, discard the Igloo Wound Irrigator appropriately.

Disposal

Igloo Wound Irrigation Shield should be disposed of according to current local, regional, state, and/or federal laws and regulations pertaining to the safe disposal of medical devices and accessories.

WARNINGS:

- Do not use if the package is open or damaged.
- For deep or puncture wounds, animal bites or serious burns, consult a healthcare professional.
- Single-use only.
- To reduce aerosolization, keep the Igloo Wound Irrigation Shield in contact with the skin during the irrigation procedure.

Medical Device Reporting: Notice to Users and/or Patients in EU: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Notice for Users and/or Patients in Switzerland: Any serious incident that has occurred in relation to the device should be reported without delay to Bionix at bionix@bionix.com and to Swissmedic (www.swissmedic.ch).

Contraindications: The Igloo Wound Irrigation Shield is not indicated for the cleansing of deep or puncture wounds.

Residual Risk: Risk associated with the use of this product has been reduced as far as possible, but the product cannot completely eliminate potential patient or user harm arising from the following:

- Harm from mechanical hazards
- Harm from misuse, or use error
- Harm from unanticipated origins

	Medical Device		Do Not Use if Package is Damaged
	Manufacturer		Do Not Reuse Single-Use Device
	Date of manufacture		Consult Instructions For Use (IFU)
	Lot Code		Caution
	Reorder Number		Warning
	Sterilized using Irradiation		Do Not Resterilize
	Prescription Only or "For Use by or on the order of a licensed medical professional"		

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