

Lighted Articulating Ear Curette™



**Precision illumination.
Behind-the-obstruction access.
Confident care.**



Why Clinicians Choose Bionix®

A smarter way to remove cerumen and foreign bodies

More than a tool, it's a purpose-built solution designed to support procedural clarity, safety and clinical efficiency.

1. Designed for safer, more precise procedures

- Integrated light and magnification give clinicians a clear view deep into the ear canal.
- The 90° articulating VersaScoop® tip allows access behind the occlusion, enabling controlled removal under direct visualization and supporting safer cerumen and foreign body removal.

2. Enhanced patient comfort and confidence

- Working behind the obstruction means less manipulation and less risk of canal irritation – improving patient experience, especially in pediatric or sensitive cases.

3. Simple, intuitive operation

- A familiar curette form factor combined with integrated illumination and magnification makes it easy for clinicians to adopt and use effectively across care settings.

4. A smart investment

- Single-use tips support infection control. The Lighted Articulating Ear Curette is also available now for use with the Rechargeable Light Source for long-term saving and value.



**Innovative Solutions
for Everyday Problems**

1670 Indian Wood Circle
Maumee, OH 43537
United States

800.551.7096 (+1 419.727.8421)

800.455.5678 (+1 419.727.8426)

bionix.com

bionix@bionix.com (international@bionix.com)

Lighted Articulating Ear Curette™



Why Behind-the-Obstruction Access Matters

- Obstructions often cannot be safely removed from the front.
- The articulating VersaScoop® tip lets clinicians get behind the material eliminating blind sweeping motions and providing direct visualization throughout the procedure.

Product Offering:

#2511 – Lighted Articulating VersaScoop (3mm tip)

Includes:

- 25 Lighted Articulating VersaScoop Tips
- 1 Disposable Light Source
- 1 Magnification Lens

#2512 – Lighted Articulating VersaScoop (3mm tip) Refill Pack

Includes:

- 25 Lighted Articulating VersaScoop Tips

#2200 – Extra Disposable Light Source

#3600 – Rechargeable Light Source

#2265 – Pack of 5 Magnification Lenses

Cerumen or Foreign Body CPT Codes:

69210 – Removal of impacted cerumen

69200 – Removal of foreign body from ear canal

30300 – Removal of foreign body from the nose

Backed by Bionix®

With more than 40 years of innovation across ear, nose and throat care, Bionix is trusted by clinicians for tools that enhance visibility, precision and patient-focused performance.



Contact Us

Contact your Bionix sales rep, distributor of choice or call 1.800.551.7096 (US)

International: Contact your distributor or call +1 419.727.8421

Lighted Articulating Curettes are Single-Use Only

Rx ONLY

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

US Patent #8,728,107

©2026 Bionix LLC | MKT-0033 | REV. 02



1670 Indian Wood Circle
Maumee, OH 43537
United States

800.551.7096 (+1 419.727.8421)
800.455.5678 (+1 419.727.8426)
bionix.com
bionix@bionix.com (international@bionix.com)