

Lighted Ear Curette™ with Magnification



**Precision illumination.
Clear visualization. Confident care.**



The Bionix® Lighted Ear Curette™ provides integrated illumination and magnification to enhance visualization of the ear canal during cerumen and foreign body removal. The design supports standard curettage technique while allowing clinicians to maintain a clear view of the treatment area.

Why Clinicians Choose Bionix

- **Integrated illumination and magnification**
Light and magnification at the curette tip support visualization of the ear canal during curettage.
- **Supports procedural control**
Integrated illumination helps clinicians maintain a consistent line of sight during routine cerumen removal.
- **Familiar technique**
Supports standard curettage techniques without requiring changes to workflow or handling.
- **Single-use curettes**
Disposable curettes support infection control protocols and consistent performance from patient to patient.

A trusted solution used in everyday practice

Bionix Lighted Ear Curettes have been used by healthcare providers for decades and are supported by a broad portfolio of tip styles designed to accommodate varying patient anatomies and cerumen presentations.

From pediatric to adult patients, clinicians can select the tip shape and size that best fits their clinical preference and the task at hand.

“Bionix Lighted Ear Curettes have really changed our practice. We can quickly and painlessly clean out ears. We love them!”

- Nancy Van Keuls, MD

Lighted Ear Curette™ with Magnification



Product offerings to match your workflow

A full range of Lighted Ear Curette tip styles is available. Standard Packs and Clinic Packs provide flexible options based on patient volume and practice needs.

Product Offering:

Includes 50 curettes per box + 1 Light Source + 1 Magnification Lens



#2210 – Lighted AngleLoop®
Angled 4mm diameter loop, similar to AngleLoop Red Safe Ear Curette™



#2240 – Lighted FlexLoop®
4mm wide tip, similar to FlexLoop White Safe Ear Curette



#2270 – Lighted VersaLoop®
3mm wide tip, similar to VersaLoop Purple Safe Ear Curette



#2220 – Lighted MicroLoop®
3mm diameter loop, similar to MicroLoop Green Safe Ear Curette



#2250 – Lighted CeraSpoon®
4mm wide tip, similar to CeraSpoon Yellow Safe Ear Curette



#2280 – Lighted CeraPik®
2.5mm concave tip bent to 45°. Great for hard, impacted cerumen and foreign bodies.



#2230 – Lighted WaveCurette™
2.5mm wave-like wide tip.



#2260 – Lighted InfantScoop®
2mm wide tip, similar to InfantScoop Blue Safe Ear Curette



#2245 – Variety Pack
Includes all tip styles except CeraPik®

Clinic Packs:

Includes 200 curettes per box + 4 Light Sources + 4 Magnification Lenses

#2201 – Lighted AngleLoop® Clinic Pack

#2202 – Lighted MicroLoop® Clinic Pack

#2203 – Lighted WaveCurette™ Clinic Pack

#2204 – Lighted FlexLoop® Clinic Pack

#2205 – Lighted CeraSpoon® Clinic Pack

#2206 – Lighted InfantScoop® Clinic Pack

#2207 – Lighted VersaLoop® Clinic Pack

#2208 – Lighted CeraPik® Clinic Pack

#2209 – Variety Clinic Pack (Includes all tip styles except CeraPik®)

Potential Reimbursement (U.S. only):

CPT Code 69210 – For potential reimbursement of impacted cerumen removal.

“The variety of sizes and tip shapes makes ear wax removal easier than ever before.”

- Carl A. Moritz, Jr., MD

Contact Us

Contact your Bionix sales rep, distributor of choice or call 1.800.551.7096 (US)


International: Contact your distributor or call +1 419.727.8421







Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

©2026 Bionix LLC | MKT-0034 | REV. 02

 **Innovative Solutions
for Everyday Problems**

 1670 Indian Wood Circle
Maumee, OH 43537
United States

 800.551.7096 (+1 419.727.8421)
 800.455.5678 (+1 419.727.8426)
 bionix.com
 bionix@bionix.com (international@bionix.com)