

Lighted Placement Tool

Bionix® Placement Tool combines illumination, magnification, and 5 depth guides to help Audiologists take less time to safely and accurately place a foam block deep in the ear.



PROBLEM

Placing a foam block deep in the ear canal is difficult without proper visualization.

SOLUTION

The Lighted Placement Tool combines the award winning technology of the Lighted Ear Curette™ with Magnification with a unique, bayonet design for an enhanced site line.

RESULTS

An elegantly slim, brilliantly lit and simply magnified device to promote proper visualization of the ear canal for more accurate foam block placement.

PRODUCT OFFERING

- #2850** - Lighted Placement Tool
(1 Light Source, 1 Magnification Lens, 15 Lighted Placement Tips)
- #2200** - Extra Light Source
- #2265** - Pack of 5 Magnification Lenses
- #2299** - Light Source Lanyard

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

800.551.7096 • www.bionix.com

©2021 Bionix LLC. | US Patent # D771815 | MKT-0036 | Rev. C