



MADE IN THE USA

Eliminate the mess of cerumen removal during irrigation

OtoClear® Ear Irrigation Tip

Safe, soft, comfortable 3-streamed-tip that prevents the risk of over-insertion



FEATURES

- Made with soft, flexible material for patient comfort
- Uses three divergent streams of water to safely loosen cerumen
- Directs water toward the ear canal wall instead of at the tympanic membrane
- Flared design prevents backsplash

BENEFITS

- Compatible with a variety of devices
- Improves efficiency and patient safety
- Reduces patient discomfort
- Reduces messy clean up of cerumen removal during irrigation



Bionix OtoClear Ear Irrigation Tip is a soft, flared, 3-streamed-tip that improves patient safety and eliminates the mess of cerumen removal during irrigation.

It's compatible with a variety of delivery devices, meeting personal preferences of healthcare professionals.

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.



OTOCLEAR TIPS ARE
SINGLE USE ONLY

Rx ONLY



OTOCLEAR TIPS ARE
MADE IN THE USA

Contact your Bionix Sales Rep, Distributor of Choice or call 1.800.551.7096 (US)
International: Contact your Distributor or call +1 419.727.8421

INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

1670 Indian Wood Circle | Maumee, OH 43537
United States | P: 800.551.7096 (+1 419.727.8421)
F: 800.455.5678 (+1 419.727.8426) | www.bionix.com
Email: bionix@bionix.com (international@bionix.com)

©2022 Bionix LLC
MKT-0041 | REV. 03

OtoClear® Ear Irrigation Tip

Product Offering



OtoClear Ear Irrigation System*

Recommended for high volume facilities

- A completely portable irrigation delivery device that can be used for procedures anywhere in or out of the office.
- Provides effective, pulsed irrigation.
- Higher volume than manual devices.
- Electronic pressure-control system with three settings.
- Proven safety and effectiveness.
(Children's Hospital of Pittsburgh Study)
- Rechargeable battery.
- 278 ml self-contained water reservoir.

#7280 - (1) Portable Electronic Irrigator,
(40) OtoClear Tips, (1) Ear Basin,
(1) 110v Charging Adapter

#7245 - (1) Portable Electronic Irrigator,
(5) OtoClear Tips, (1) Ear Basin,
(1) 110v Charging Adapter

#7282 - (1) Portable Electronic Irrigator,
(40) OtoClear Tips, (1) Ear Basin,
(1) 220v Charging Adapter

#7247 - (1) Portable Electronic Irrigator,
(5) OtoClear Tips, (1) Ear Basin,
(1) 220v Charging Adapter

***PRODUCT AVAILABILITY VARIES BY COUNTRY**



OtoClear SprayWash Kit

- Manual, intermittent lavage.
- Delivers greater volume compared to syringes or other manual devices.
- Efficient design allows for two-handed procedure.
- 500ml bottle volume.

#7290 - (1) SprayWash Bottle,
(20) OtoClear Tips, (1) Ear Basin

#7295 - (1) SprayWash Bottle,
(5) OtoClear Tips, (1) Ear Basin



OtoClear Ear Irrigation Tips

- Directs water to the ear canal wall, with three divergent streams instead of directly at the tympanic membrane.
- Flared design prevents backslash and over-insertion.

#7200 - OtoClear Ear Irrigation Tips
(40 Tips per Box)

#7202 - OtoClear Ear Irrigation Tips
(20 Tips per Box)



OtoClear Aquabot® Kit

- Delivers greater volume compared to syringes or other manual devices.
- Efficient design allows for two-handed procedure.
- Provides continuous stream without manually pumping.
- 650ml bottle volume.

#7270 - (1) Aquabot Bottle,
(20) OtoClear Tips, (1) Ear Basin

#7275 - (1) Aquabot Bottle,
(5) OtoClear Tips, (1) Ear Basin

Aquabot is registered under Lunatech®, used under permission of Lunatech, Inc. All rights reserved.
The Portable Electronic Irrigator and Aquabot Bottle are not manufactured in the USA.

Contact your Bionix Sales Rep, Distributor of Choice or call 1.800.551.7096 (US)
International: Contact your Distributor or call +1 419.727.8421