

OtoClear® Ear Irrigation Delivery Devices



A OtoClear Ear Irrigation System

B OtoClear SprayWash Kit

C OtoClear Aquabot® Kit

Description	Offering	Device System	Pressure Control Settings	Water Reservoir	Recommended Facilities
A OtoClear Ear Irrigation System PRODUCT AVAILABILITY VARIES BY COUNTRY	#7280 - (1) Portable Electronic Irrigator, (1) Ear Basin, (40) OtoClear Tips, (1) 110v Charging Adapter #7245 - (1) Portable Electronic Irrigator, (1) Ear Basin, (5) OtoClear Tips, (1) 110v Charging Adapter #7282 - (1) Portable Electronic Irrigator, (1) Ear Basin, (40) OtoClear Tips, (1) 220v Charging Adapter #7247 - (1) Portable Electronic Irrigator, (1) Ear Basin, (5) OtoClear Tips, (1) 220v Charging Adapter	Automatic	Low Medium High	278ml	High Volume Facilities
B OtoClear SprayWash Kit	#7290 - 1 SprayWash Bottle, 20 OtoClear Tips, 1 Ear Basin #7295 - 1 SprayWash Bottle, 5 OtoClear Tips, 1 Ear Basin	Manual	Volume Based System	500ml	Low Volume Facilities
C OtoClear Aquabot Kit	#7270 - 1 Aquabot Bottle, 20 OtoClear Tips, 1 Ear Basin #7275 - 1 Aquabot Bottle, 5 OtoClear Tips, 1 Ear Basin	Manual	Stealth Mode Pulse Mode	650ml	Low Volume Facilities

Contact your Bionix Sales Rep, Distributor of Choice or call 1.800.551.7096 (US)
International: Contact your Distributor or call +1 419.727.8421

©2022 Bionix LLC
MKT-0076 | REV. 03

INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

1670 Indian Wood Circle | Maumee, OH 43537
United States | P: 800.551.7096 (+1 419.727.8421)
F: 800.455.5678 (+1 419.727.8426) | www.bionix.com
Email: bionix@bionix.com (international@bionix.com)