

OtoClear® Portable Electronic Irrigator

A versatile solution for an automated ear irrigation device.

The Bionix OtoClear Portable Electronic Irrigator is designed for safe and effective cerumen removal when used with our patented OtoClear Tips. The Tip's unique, flared design prevents over-insertion while prohibiting any backsplash. The three-stream water flow creates a circular clearing method around the ear canal wall, directing water away from the tympanic membrane. The result is a safe, simple and comfortable approach to ear wax removal.

PRODUCT OFFERING



Recommended for high volume facilities

Bionix OtoClear Portable Electronic Irrigator Features:

- Rechargeable
- Fully portable
- Adjustable flow settings
- Pulsed irrigation
- Electronic pressure control system

OtoClear Ear Irrigation System Offerings:

#7280 - (1) Portable Electronic Irrigator,
(1) Ear Basin, (40) OtoClear Tips,
(1) 110v Charging Adapter

#7245 - (1) Portable Electronic Irrigator,
(1) Ear Basin, (5) OtoClear Tips,
(1) 110v Charging Adapter

#7282 - (1) Portable Electronic Irrigator,
(1) Ear Basin, (40) OtoClear Tips,
(1) 220v Charging Adapter

#7247 - (1) Portable Electronic Irrigator,
(1) Ear Basin, (5) OtoClear Tips,
(1) 220v Charging Adapter

#7200 - Box of 40 OtoClear Ear Irrigation Tips

#7202 - Box of 20 OtoClear Ear Irrigation Tips



PRODUCT AVAILABILITY VARIES BY COUNTRY

Contact your Bionix Sales Rep, Distributor of Choice or call 1.800.551.7096 (US)
International: Contact your Distributor or call +1 419.727.8421

©2022 Bionix LLC
MKT-0083 | REV. 03

INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

1670 Indian Wood Circle | Maumee, OH 43537
United States | P: 800.551.7096 (+1 419.727.8421)
F: 800.455.5678 (+1 419.727.8426) | www.bionix.com
Email: bionix@bionix.com (international@bionix.com)