



KIT NOW AVAILABLE!



Game Changer for Suture Removal

The only disposable tool that safely cuts sutures from outside the suture loop with a guarded-blade. The **SutureOut** tool maintains the integrity of the newly healed wound while minimizing patient anxiety.

Adhesive Pad



Guarded-blade

✓ Top Reasons to Use The **SutureOut** Kit

- Fast, simple suture removal
- Better control of suture removal tools and process
- No exposed blade
- Sterile, disposable and environmentally friendly
- Unique adhesive pad for collecting and counting removed sutures
- Keeps wound integrity, even with tight sutures
- Less tugging, height and tension needed to cut the suture
- Designed with an outside-in approach
- No reprocessing of metal tools



The **SutureOut** tool solves clinical problems when faced with thin skin, tightly tied knots, embedded stitches, or patients experiencing fear and anxiety.

The **SutureOut** tool has a minimalistic design to reduce cost and medical waste associated with the use of scissors and standard metal suture removal kits.



**Innovative Solutions
for Everyday Problems**

1670 Indian Wood Circle
Maumee, OH 43537
United States

800.551.7096 (+1 419.727.8421)

800.455.5678 (+1 419.727.8426)

bionix.com

bionix@bionix.com (international@bionix.com)

Product Offering:

KIT NOW AVAILABLE!

#8520 – SutureOut Suture Removal Kit (50 Kits per Case)

Each Kit contains:

- (1) SutureOut suture removal tool
- (1) 4" x 4" gauze sponge
- (1) Stainless steel Adson Forceps without teeth

#8500 – SutureOut Suture Removal Tool (50 Tools per Box)



Suture Removal CPT Codes:

15851 Removal of sutures or staples requiring anesthesia (e.g., General anesthesia, moderate sedation)

15853 Removal of sutures or staples not requiring anesthesia (List separately in addition to E/M code)

15854 Removal of sutures and staples not requiring anesthesia (List separately in addition to E/M code)

“ I previously thought using curved scissors was best practice for suture removal. After using **SutureOut**, I realized it is a much better option. I no longer have to pull up on the suture, which can cause pain for my patient. ”

Pediatrician, Perrysburg, Ohio

Contact Us

Contact your Bionix sales rep, distributor of choice
or call 1.800.551.7096 (US)

International: Contact your distributor or call
+1 419.727.8421

 **Rx ONLY**  **STERILE**  **R**

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

©2025 Bionix LLC | Patent #US9532779B2 | MKT-0296 | REV. 02



**Innovative Solutions
for Everyday Problems**

1670 Indian Wood Circle
Maumee, OH 43537
United States

800.551.7096 (+1 419.727.8421)
800.455.5678 (+1 419.727.8426)
bionix.com
bionix@bionix.com (international@bionix.com)