

Rechargeable Light Source

for ClearLook® Lighted Ear Curettes™



**Precision, Efficiency and
Sustainability for Safer Ear Care**



✓ The Rechargeable Light Source for ClearLook Lighted Ear Curette offers:

- **Trusted Quality:** Innovative tools from Bionix® – the leader in ear, nasal and throat care products.
- **Rechargeable Convenience:** Includes a USB-C charging cord; an optional charging base can be purchased.
- **Cost-Effective:** Affordable solution with long-term curette refill savings and improved procedure efficiency.
- **Eco-Friendly:** Rechargeable design reduces waste and costs.
- **No Blind Curettage:** Light and magnification provide a clear view of the ear canal.
- **Streamlined Workflow:** No need to switch between otoscope and curette.
- **Precision and Safety:** Enhanced visibility reduces risk of injury and unseen complications
- **Patient Comfort:** Faster, more precise and less intrusive procedures.
- **Seamless Integration:** Compatible with preferred curette tips, intuitive to use.

Premium Performance in Your Hands


The Bionix® Rechargeable Light Source for ClearLook® Lighted Ear Curette™ enhances cerumen removal with improved visibility, precision and eco-friendly technology. Its rechargeable design reduces waste, while compatibility with familiar and preferred curette tip styles ensures a seamless, cost-effective upgrade for clinicians.

Rechargeable Light Source

for ClearLook® Lighted Ear Curettes™



Product Offering:



#3645 – Rechargeable Light Source with ClearLook Variety Pack - Starter System


- (1) Rechargeable Light Source
- (50) Assorted ClearLook Lighted Ear Curettes (All 6 Styles)
- (1) Magnification Lens
- (1) USB-A to USB-C charging cable
- (1) Charging block



#3600 – Rechargeable Light Source

- (1) Rechargeable Light Source
- (1) Magnification Lens
- (1) USB-A to USB-C charging cable
- (1) Charging block

Also compatible with:



#3410 – ClearLook® AngleLoop® Refill Pack (50 Curettes) Angled 4mm diameter loop, similar to AngleLoop Red Safe Ear Curette™




#3420 – ClearLook® MicroLoop® Refill Pack (50 Curettes) 3mm diameter loop, similar to MicroLoop Green Safe Ear Curette



#3430 – ClearLook® VersaLoop® Refill Pack (50 Curettes) 3mm wide tip, similar to VersaLoop Purple Safe Ear Curette



#3440 – ClearLook® FlexLoop® Refill Pack (50 Curettes) 4mm wide tip, similar to FlexLoop White Safe Ear Curette



#3450 – ClearLook® CeraSpoon® Refill Pack (50 Curettes) 4mm wide tip, similar to CeraSpoon Yellow Safe Ear Curette



#3460 – ClearLook® InfantScoop® Refill Pack (50 Curettes) 2mm wide tip, similar to InfantScoop Blue Safe Ear Curette



#3445 – ClearLook® Variety Pack Refill Pack (50 Curettes) An assortment of all 6 styles of ClearLook Lighted Ear Curettes



#3695 – Multi-Port Charging Base
1 charging base with USB-C to USB-C cable
(charges up to 4 Rechargeable Light Sources; not included with base.)

Part of the **OneLight™**
Illumination Platform

Contact Us

Contact your Bionix sales rep, distributor of choice or call 1.800.551.7096 (US)

International: Contact your distributor or call +1 419.727.8421

Rx ONLY



ClearLook Lighted Ear Curettes are Single Use Only

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

ClearLook® Lighted Ear Curette™ US Patent #s: 12402974, D1036670, D1090838

©2026 Bionix LLC | MKT-0314 | REV. 02



**Innovative Solutions
for Everyday Problems**

1670 Indian Wood Circle
Maumee, OH 43537
United States

800.551.7096 (+1 419.727.8421)

800.455.5678 (+1 419.727.8426)

bionix.com

bionix@bionix.com (international@bionix.com)