

# Lighted MicroSuction System



NEW



**Precision illumination.**  
**Quiet operation. Confident care.**

## Why Clinicians Choose Bionix

A better system to remove cerumen and foreign bodies

More than a tool, it's a smarter system built to support procedural clarity, patient comfort and clinical efficiency.

### 1. Designed for safer, more precise procedures

Integrated light and magnification give clinicians a clear view of the ear or nasal passage, improving visibility, supporting accuracy and reducing the chance of complications.

### 2. Enhanced patient comfort

The hospital-grade suction pump, developed by a pediatric ENT, operates at 54 dB and offers three levels of effective suction so you can match the power to the procedure.

### 3. Compact and easy to use pump

The cordless pump is rechargeable, portable and lightweight. Its small footprint makes it ideal for use in all patient care settings.

### 4. A smart investment

A complete system that supports clinical performance and adds value through streamlined efficiency.

### 5. Versatile for ear and nasal procedures

Supports both cerumen removal and foreign body extraction.




HOW TO  
USE VIDEO



## Why Quiet Matters

Noise levels impact both patients and providers.



54 dB = Bionix MicroSuction System Pump – equivalent to a soft conversation

69 dB = standard ENT vacuum pump

A 15 dB reduction makes the procedure more comfortable for patients and less stressful for your team.



Innovative Solutions  
for Everyday Problems

1670 Indian Wood Circle  
Maumee, OH 43537  
United States

800.551.7096 (+1 419.727.8421)

800.455.5678 (+1 419.727.8426)

bionix.com

bionix@bionix.com (international@bionix.com)

# Lighted MicroSuction System



## The Bionix Lighted MicroSuction System

Combining bright illumination, magnification and effective suction in one compact, cordless solution.

Designed to improve visibility and streamline procedures, this complete system supports confident removal of cerumen and foreign bodies while enhancing the overall clinical experience.

### Product Offering:

#### #2650 - Lighted MicroSuction System

- (1) Rechargeable, Hospital-grade Suction Pump
- (10) Lighted Suction Handles
- (1) Rechargeable Light Source
- (1) Magnification Lens
- (1) Flexible Suction Tubing (39.4")
- (1) Cerumen Catcher
- (2) Charging Cables
- (2) Charging Blocks

#### #2651 - (25) Lighted Suction Handle Refill Pack

#### #2652 - (1) Replacement Cerumen Catcher

#### #2653 - (1) Replacement Suction Tubing

### Cerumen or Foreign Body CPT Codes:

69210 – Removal impacted cerumen (separate procedure), 1 or both ears.

69200 – Removal foreign body from external auditory canal; without general anesthesia.

30300 – Removal foreign body, intranasal; office type procedure.



### Backed by Bionix

With more than 40 years as the leader of ear, nose and throat innovation, Bionix is trusted by clinicians across the continuum of care. The Lighted MicroSuction System delivers the visibility, precision and patient-focused design that providers value most.

## Contact Us

Contact your Bionix sales rep, distributor of choice or call 1.800.551.7096 (US)

International: Contact your distributor or call  
+1 419.727.8421

### Lighted Suction Handles are Single-Use Only

#### ONLY

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

Lighted Suction Handle US Patent # 9,114,202

©2025 Bionix LLC | MKT-0321 | REV. 01



Innovative Solutions  
for Everyday Problems

1670 Indian Wood Circle  
Maumee, OH 43537  
United States

800.551.7096 (+1 419.727.8421)  
800.455.5678 (+1 419.727.8426)  
bionix.com  
bionix@bionix.com (international@bionix.com)