

OtoClear® Powered by Aquabot®



Streamlined Ear Irrigation Designed for Clinical Efficiency

The Bionix® OtoClear® Aquabot® System is an advanced, easy-to-use solution designed to simplify and enhance ear irrigation in busy medical environments. This innovative system pairs the precision of the OtoClear Ear Irrigation Tip with the efficiency of the Aquabot delivery device, offering a streamlined approach to cerumen removal for healthcare professionals.

With the Aquabot System, a single touch activates a continuous and controlled stream of water, ensuring a comfortable and thorough cerumen removal process for patients. The OtoClear Tip's three-stream technology, directs water away from the eardrum, promoting patient safety. Exit portals on the OtoClear Tip help minimize splash-back.

Ideal for physician offices, urgent care centers, and healthcare facilities, the OtoClear Aquabot System improves workflow efficiency while delivering reliable, repeatable results. Designed for ease of use and patient comfort, this system enables healthcare providers to focus on what matters most – providing exceptional care.

Transform your ear irrigation process with the Bionix OtoClear Aquabot System: a safe and more efficient way to meet the demands of your practice.

Designed to Simplify Ear Irrigation in Busy Clinical Settings

Product Offering:




OtoClear® Aquabot® System

- Delivers greater volume compared to syringes or other manual devices.
- Efficient design allows for two-handed procedure.
- Provides continuous stream without manually pumping.
- 650ml bottle volume.

#7270 – (1) Aquabot Bottle, (20) OtoClear Tips, (1) Ear Basin

#7275 – (1) Aquabot Bottle, (5) OtoClear Tips, (1) Ear Basin



OtoClear Ear Irrigation Tips

- Directs water toward the ear canal walls through three divergent streams rather than directly at the tympanic membrane.
- Flared design prevents over-insertion.

#7200 – OtoClear Ear Irrigation Tips (40)

#7202 – OtoClear Ear Irrigation Tips (20)

Backed by Bionix

With more than 40 years of experience developing clinician-trusted medical devices, Bionix delivers practical solutions designed to support everyday patient care. The OtoClear Aquabot Ear Irrigation System reflects Bionix's commitment to performance, reliability and clinician confidence across care settings.

Contact Us

Contact your Bionix sales rep, distributor of choice
or call 1.800.551.7096 (US)

International: Contact your distributor or call
+1 419.727.8421



OTOCLEAR TIPS ARE
SINGLE USE ONLY

Rx ONLY



Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

Aquabot is registered under Lunatech®, used under permission of Lunatech, Inc. All rights reserved.

©2026 Bionix LLC | MKT-0326 | REV. 01



1670 Indian Wood Circle
Maumee, OH 43537
United States

800.551.7096 (+1 419.727.8421)
800.455.5678 (+1 419.727.8426)
bionix.com
bionix@bionix.com (international@bionix.com)