ClearLook® Lighted Ear Curettes™ with OneLight™ Rechargeable Light Source





See where you're working. Perform confidently.

Designed for visibility and control

ClearLook® Lighted Ear Curettes™ bring illumination, magnification and thoughtful instrument design together to support controlled cerumen and foreign body removal from the ear canal.

A natural, pencil-like grip keeps finger placement outside the line of sight, helping clinicians maintain their view while working. Six loop and scoop tip styles provide options for instrument selection based on the patient and clinical assessment.

Enhanced by the OneLight™ Rechargeable Light Source

Pair ClearLook with the OneLight™ Rechargeable Light Source to bring focused illumination to the curette during use. The Magnification Lens complements the light to support your view of the ear canal.

With up to two hours of continuous use per charge, the Rechargeable Light Source is designed for continued use across procedures. Retain and clean the light source according to its instructions, then replenish your single-use ClearLook curettes as needed.

Backed by Bionix Innovation

ClearLook reflects more than 40 years of Bionix innovation in ear, nose and throat care, bringing clinician-focused design and OneLight illumination together in the exam room.


Why ClearLook Matters

- **Visibility at the point of use:** Illumination and magnification support visualization of the ear canal.
- **Purposeful handle design:** Secure finger placement outside the line of sight and a natural, pencil-like grip support instrument control.
- **Six tip styles:** Select from AngleLoop®, MicroLoop®, VersaLoop®, FlexLoop®, CeraSpoon® and InfantScoop®.
- **Single-use curettes:** Use a new ClearLook curette for each patient and discard after use.

Choose your ClearLook starter system or refill packs. See Page 2 for tip configurations and ordering information.

ClearLook® Lighted Ear Curettes™ with OneLight™ Rechargeable Light Source

Product Offering:



#3645 – OneLight™ Rechargeable Light Source with ClearLook® Variety Pack – Starter System
- (1) OneLight™ Rechargeable Light Source
- (50) Assorted ClearLook Lighted Ear Curettes™ (All 6 Styles)
- (1) Magnification Lens
- (1) USB-A to USB-C charging cable
- (1) Charging block



#3600 – OneLight™ Rechargeable Light Source
- (1) OneLight™ Rechargeable Light Source
- (1) Magnification Lens
- (1) USB-A to USB-C charging cable
- (1) Charging block

Also compatible with:



#3410 – ClearLook® AngleLoop® Refill Pack (50 Curettes) Angled 4mm diameter loop, similar to AngleLoop Red Safe Ear Curette™



#3420 – ClearLook® MicroLoop® Refill Pack (50 Curettes) 3mm diameter loop, similar to MicroLoop Green Safe Ear Curette



#3430 – ClearLook® VersaLoop® Refill Pack (50 Curettes) 3mm wide tip, similar to VersaLoop Purple Safe Ear Curette



#3440 – ClearLook® FlexLoop® Refill Pack (50 Curettes) 4mm wide tip, similar to FlexLoop White Safe Ear Curette



#3450 – ClearLook® CeraSpoon® Refill Pack (50 Curettes) 4mm wide tip, similar to CeraSpoon Yellow Safe Ear Curette



#3460 – ClearLook® InfantScoop® Refill Pack (50 Curettes) 2mm wide tip, similar to InfantScoop Blue Safe Ear Curette



#3445 – ClearLook® Variety Pack Refill Pack (50 Curettes) An assortment of all 6 styles of ClearLook Lighted Ear Curettes



#3695 – Multi-Port Charging Base
1 charging base with USB-C to USB-C cable
(charges up to 4 Rechargeable Light Sources: not included with base.)

▶ Explore Bionix Clinical Education

Contact Us

Contact your Bionix sales rep, distributor of choice or call 1.800.551.7096 (US)

International: Contact your distributor or call +1 419.727.8421

℞ ONLY

ClearLook Lighted Ear Curettes are Single Use Only

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions.

ClearLook® Lighted Ear Curette™ US Patent #s: 12402974, D1036670, D1090838

©2026 Bionix LLC | MKT-0314 | REV. 03

bionix®
Innovation With Impact

1670 Indian Wood Circle
Maumee, OH 43537
United States

800.551.7096 (+1 419.727.8421)
800.455.5678 (+1 419.727.8426)
bionix.com
bionix@bionix.com (international@bionix.com)